

4 WHEEL ROLLATORS (£15*) £18 PER WEEK TO HIRE

STANDARD ROLLATOR (*120) £144.00

Lightweight Aluminium, less than 7kg
Adjustable height handles 32"-37"
Adjustable height seat 19"-23"
Arthritic friendly, looped cable brakes
Anatomical handgrips
Padded flip-up seat
Underseat storage bag
Available in blue or red
Max user weight 21st / 135kg



NITRO SL ROLLATOR (*155) £186.00

Ultra Lightweight Aluminium 5.5kg
Adjustable height handles 32.7-36.6"
Adjustable height backrest
Arthritic friendly, looped cable brakes
Built in stick holder
Removable storage bag
Clip holds frame closed when folded
Available in blue or grey
Max user weight 136kg/21.4st

SATURN ROLLATOR (*199.99) £239.99

Lightweight & modern 7.3kg
8" front soft swivel wheels
Adjustable height handles 34"-41"
Arthritic friendly, looped cable brakes
Built in stick & drink holder
Removable storage bag
Available in black or champagne
Max user weight 22st/136kg



4 WHEEL ROLLATORS (£15*) £18 PER WEEK TO HIRE

TAiMA M-GT (*261.66) £313.99

Ultra lightweight 5.5kg
 As light as a 3-wheeled walker
 Adjustable height handles
 Soft ergonomic handles
 Walking stick holder included
 Large storage bag
 Arthritic friendly, looped cable brakes
 Front & side reflectors
 Small & Medium models available
 Max user weight 23st/150kg



TAiMA XC (*429.99) £515.99

Pneumatic tyres & XL wheels
 Brilliant for multi surface outdoor use
 Adjustable handle height
 Seat height 61cm
 Total weight 9.14kg
 Max user weight 23st/150kg

TOPRO TROJA 5G (*356.00) £427.20

Back rest included
 Adjustable height handles
 Available in 8 modern colours
 Easy folding, lightweight 7.8kg
 Soft grip / ergo grip handles
 Huge range of accessories
 7 year warranty
 2 types of wheels available
 2 frame sizes available
 Max user weight 23.5st / 150kg



HIRE ONLY: Local BN16 & BN17: Delivery & collection charge = (*20.00) £24.00
HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (*40.00) £48.00
HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (*60.00) £72.00

***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**

Try before you buy:
If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.

4 WHEEL ROLLATORS (£15*) £18 PER WEEK TO HIRE

EZ FOLD-N-GO ROLLATOR (*249.99) £299.99

Folds up to 3x smaller than average
Portable & lightweight 6.4kg
Adjustable height handles
Padded handles & locking brakes
2 colours available
Underseat storage bag
Organiser for on the go storage
Max user weight 17.5st / 113kg



EZ Fold-N-Go Walker (*114.99) £137.99

Super lightweight 3.6kg
2x wheels & 2x easy-glide feet
No seat & no brakes
Single handed folding
Compact & portable
Storage pouch
Adjustable hand height 81-97cm
Available in blue or walnut
Max user weight 28st / 180kg

© HANDBELT 2022

"LETS GO INDOOR" ROLLATOR (*264.99) £317.99

Lightweight 6.5kg
Adjustable height single handle
Very narrow: 54.5cm wide
Full length soft handlebar
Full length one-handed brake
Removable tray
Tray colour may vary
Undertray storage bag
Foldable



***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**

Try before you buy:

If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.

4 WHEEL ROLLATORS (£15*) £18 PER WEEK TO HIRE



HANDY TROLLEY (*89.00) £106.80

Height adjustable handles
Locking brakes on handles
2x Removable trays
Non-folding
Indoor use only

RUTLAND TROLLEY (*72.50) £87.00

Adjustable height 84-99cm
Two trays
Non-folding
No brakes
Perfect for moving items around the house
Max user weight 28st /180kg



AirGo Fusion 2-in-1 (*205) £246.00

Rollator and transport chair
Easily converts between 2 uses
Swing-around leg & back rests
Adjustable height 78.5-92.5cm
Padded seat & backrest
Removable storage bag
67cm width
Weight 9.5kg
Max user weight 21st/ 136kg

HIRE ONLY: Local BN16 & BN17: Delivery & collection charge = (*20.00) £24.00
HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (*40.00) £48.00
HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (*60.00) £72.00

***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**

Try before you buy:
If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.

PREMIUM 4 WHEEL ROLLATORS - Not for hire, but can be trialed in the showroom.

ROLLZ FLEX 2 (*499.00) £598.80

Stylish, designer rollator, 7.7kg
Large, water resistant, shopping bag
Adjustable angle & height handlebar
Less vibrations from foam-filled tyres
Small turning circle
Leather look seat
Available in 4 colours.
Range of accessories available
Max user weight: 19.5st / 125kg



AEROWALK CARBON ROLLATOR (340) £480

Lightest rollator available:

4.6kg / 10.1lbs without accessories

5.1kg / 11.2lbs with accessories

Premium quality carbon fibre frame 23.6" width

Modern, stylish design with wide back strap

Easy adjustment handle height 31.5"-37.4"

Ergonomic handles. Seat Height 22.4".

Compact folding with lock to hold closed

All round reflectors

Available in red or white

Max user weight: 23.6st / 150kg



***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**