

TOPAK CARE SUPPLIES LTD.

Tel: 01903 783000

www.topak.co.uk



Fax: 01903 771477

sales@topak.co.uk

UNITS 1-2 BROOKSIDE AVENUE, RUSTINGTON, WEST SUSSEX BN16 3LF



MOBILE COMMODE

Chromium plated steel
Detachable Arm Rests
4x Braked Castors

Depth - 58cm
Width - 52cm
Height - 98cm
Seat Height - 55cm / 22"
Seat Width - 44cm
Seat Depth - 46cm / 18"
Weight - 13.1kg / 28.8lbs
Max User Weight 155kg / 341lbs

To Buy =(*125.00) £150.00
To Hire =(*10.00) £12.00 per week
Commode pan & lid **MUST** be bought = (*6.20) £7.44



COMMODE/SHOWER CHAIR

Plastic coated steel
Arms tilt up
4 Braked Castor
Adjustable Footrest
52cm
64cm
Between arms – 54cm
50cm
42cm
42cm
15kg / 32lbs
160kg / 25st

(*210) £252.00
(*15.00) £18.00 per week

Zone1 (under 20 miles): Delivery & Collection Fee= (*20.00) £24.00
Zone 2 (20-40 miles) Delivery & Collection Fee= (*30.00) £36.00

This item is eligible for *VAT relief* if you suffer a chronic medical condition.

Try before you buy:

If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.

TOPAK CARE SUPPLIES LTD.

Tel: 01903 783000

www.topak.co.uk



Fax: 01903 771477

sales@topak.co.uk

UNITS 1-2 BROOKSIDE AVENUE, RUSTINGTON, WEST SUSSEX BN16 3LF

STACKABLE STATIC COMMODE



1)



2)

Stackable design for easy of storage. Completely removable upholstery for cleaning & storage. Splayed legs for added stability, supplied with removable commode pan & lid, fire retardant upholstery. Extendable legs option available (2)
Max user weight 160kg (25st)

1) Standard - To buy = (*33.00) £39.60

2) Adjustable Height - To buy = (*43.00) £51.60

To hire either = (*10.00) £12.00 per week

Commode pan & lid **MUST** be bought = (*5.20) £6.24

Zone1 (under 20 miles): Delivery & Collection Fee= (*20.00) £24.00

Zone 2 (20-40 miles) Delivery & Collection Fee= (*30.00) £36.00

This item is eligible for *VAT relief* if you suffer a chronic medical condition.

Try before you buy:

If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.