

TOPAK CARE SUPPLIES LTD.

Tel: 01903 783000

www.topak.co.uk



Fax: 01903 771477

sales@topak.co.uk

UNITS 1-2 BROOKSIDE AVENUE, RUSTINGTON, WEST SUSSEX BN16 3LF

REPOSE CONTUR

- Riser recliner pressure redistribution overlay
- Combines enhanced comfort ... with pressure area care

- New Repose Contur is a pressure redistribution cushion overlay for Riser Recliner chairs utilising the Repose tried and trusted inflatable air cell and 'smart valve' technologies.

- Offering pressure area care to the seat, back and leg sections, it secures to the chair easily and safely with straps and thus it can

either be fitted to new chairs or be a straightforward retro-fit for users who already have a chair and require a specialist support surface.



- **Product Features**

- Adjustable straps for a secure fit on a range of chair sizes.
- Repose tried and trusted single air cell and smart valve technologies.
- Helps protect the cover material on the chair.
- Provides pressure redistribution across the whole surface for the at-risk areas of the body - sacrum, buttocks, scapulae, vertebrae, shoulders, back of head.
- Sophisticated technology in a simple to use recliner chair overlay cushion.
- Secure strap system allows for straightforward retro-fit or new chair use.
- Easy to clean and maintain.
- Cost effective method to provide pressure area care to a recliner chair.

- **Specification**

- Contents: 1 x Repose Contur Inflatable Inner Section, 1 x Contur Cover and 1 x Large Pump
- Inflated dimensions: L 1780mm x W 575mm
 - Inflated Back Section dimensions: L 650mm x W 575mm
 - Inflated Seat Section dimensions: L 500mm x W 575mm
 - Inflated Leg Section dimensions: L 520mm x W 575mm
- Max patient weight 139kg
-

To Buy = £299 + VAT or To Hire = £10 + VAT per week
Chronically ill or disabled users do not have to pay VAT