

TOPAK CARE SUPPLIES LTD.

Tel: 01903 783000

www.topak.co.uk



Fax: 01903 771477

sales@topak.co.uk

UNITS 1-2 BROOKSIDE AVENUE, RUSTINGTON, WEST SUSSEX, BN16 3LF

LIGHTWEIGHT WHEELCHAIRS (£15*) £18 PER WEEK TO HIRE

DRIVE ULTRA LIGHTWEIGHT WHEELCHAIRS (£245*) £294

Ultra lightweight aluminium frame, fully folding.

17" or 20" seat width options

Models available as transit or self propelling

Handy pocket with tools on back

Adjustable height footrests.

Desk style armrests

Puncture proof wheels

Easily fits in the back of a car.

Weighs 12kg / 26lbs. 18 stone user weight limit

Quick release wheels on self propelling model. Attendant brakes only on transit models



ENIGMA XS ALUMINIUM WHEELCHAIRS (£245*) £294

Lightweight aluminium frame, 18" or 20" seat width options
Deep comfort, tension adjustable back and seat canvasses.

Comes with deluxe seat and back cushions.

Handy pocket with tools on back

Adjustable height footrests with heel straps.

Flip up, desk style armrests

Puncture proof wheels. Easily fits in the back of a car.

Height adjustable anti-tip wheels.

Adjustable rear wheel position.

Weighs 32lbs in carry format. 18 stone / 115kg weight limit.

ROMA LIGHTWEIGHT WHEELCHAIR (£195*) £234

Height adjustable footrests with heel strap

18" seat width. Loop cable brakes.

Padded, full length, fixed arms

Epoxy coated silver tone finish

30cm rear wheel & 20cm swivel front wheels

Total weight only 11kg / 24.5lbs

Maximum user weight 110kg / 17.3 stone



Zone 1 (up to 20 miles): Delivery & collection charge = (*20.00) £24.00

Zone 2 (20-40 miles): Delivery & collection charge = (*30.00) £36.00

***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**

Try before you buy:

If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.