

4 WHEEL ROLLATORS (£15*) £18 PER WEEK TO HIRE

STANDARD ROLLATOR (*120) £144.00

Lightweight Aluminium, less than 7kg
Adjustable height handles 32"-37"
Adjustable height seat 19"-23"
Arthritic friendly, looped cable brakes
Anatomical handgrips
Padded flip-up seat
Underseat storage bag
Available in blue or red



SATURN ROLLATOR (*199.99) £239.99

Lightweight & modern 7.3kg
8" front soft swivel wheels
Adjustable height handles 34"-41"
Arthritic friendly, looped cable brakes
Built in stick & drink holder
Removable storage bag
Available in black or champagne
Max user weight 22st/136kg

TAiMA M-GT (*245.83) £294.99

Ultra lightweight 5.5kg
As light as a 3-wheeled walker
Adjustable height handles
Soft ergonomic handles
Walking stick holder included
Large storage bag
Arthritic friendly, looped cable brakes
Front & side reflectors
Small & Medium models available



4 WHEEL ROLLATORS (£15*) £18 PER WEEK TO HIRE



EZ FOLD-N-GO ROLLATOR (*237.99)

£285.59

Folds up to 3x smaller than average rollators

Portable & lightweight 6.4kg

Adjustable height handles

Padded handles & locking brakes

2 colours available

Underseat storage bag

Organiser for on the go storage

“LETS DREAM” INDOOR ROLLATOR (*319.99) £383.99

Lightweight 6.5kg

Adjustable height single handle

Very narrow: 57cm wide

Full length soft handlebar

Full length one-handed brake

Removable tray

Undertray storage bag

Foldable



HANDY TROLLEY (*85.00) £102.00

Height adjustable handles

Locking brakes on handles

2x Removable trays

Non-folding

Indoor use only

Also available without brakes.



HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (*20.00) £24.00

HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (*30.00) £36.00

***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**

Try before you buy:

If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.

PREMIUM 4 WHEEL ROLLATORS

ROLLZ FLEX 2 (*499.00) £598.80

Stylish, designer rollator, 7.7kg
Large, water resistant, shopping bag
Adjustable angle & height handlebar
Less vibrations from foam-filled tyres
Small turning circle
Leather look seat
Available in 4 colours.
Range of accessories available
Max user weight: 125kg



SALJOL CARBON ROLLATOR (*589.00) £706.80

Ultra Lightweight at 6kg
Integrated drag brake – stops roll-away
All terrain, large, shock absorbing wheels
Cork, natural feel handles
Visible in the dark – over 30 reflectors
Cane holder & removable bag
Available in 3 colours
Large range of accessories available
Seat Heights: 54cm or 62cm
Max user weight: 150kg

SALJOL ALLROUND ROLLATOR (*479.00) £574.80

Very lightweight 6.6kg
Extra large all terrain wheels 25cm
Converts to single handed brake
Upgraded back rest included
Adjustable height natural cork handles
Cane holder comes as standard
Colours: Pearl Grey or Magic Black
Seat Heights: 54cm or 62cm
Max user weight 150kg



Our premium rollators are not available for hire, but you are welcome to visit our showroom to try them out.

***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**

PREMIUM 4 WHEEL ROLLATORS

LET'S GO OUT (*329.99) £395.99

Lightweight at 7kg
Large sprung wheels
Shaped rubber handles
Height adjustable handles
Generous leather seat
Removable bag & backrest
Rear splash guards
Max user weight 127kg.



TOPRO TROJA 5G (*356.00) £427.20

Back rest included
Adjustable height handles
Available in 4 modern colours/designs
Easy folding, lightweight 7.8kg
Soft grip / ergo grip handles
Huge range of accessories
7 year warranty
3 types of wheels available
2 frame sizes available
Max user weight 150kg.

LET'S MOVE (*314.99) £377.99

Ultra Lightweight at 5.4kg
Simple, elegant & sturdy.
Height adjustable handles
Soft ergonomic handles
Larger wheels to absorb bumps
Folds to just 20cm wide
Removable bag & backrest



Our premium rollators are not available for hire, but you are welcome to visit our showroom to try them out.

***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**