

4 WHEEL ROLLATORS (£15*) £18 PER WEEK TO HIRE

STANDARD ROLLATOR (*120) £144.00

Lightweight Aluminium, less than 7kg
Adjustable height handles 32"-37"
Adjustable height seat 19"-23"
Arthritic friendly, looped cable brakes
Anatomical handgrips
Padded flip-up seat
Underseat storage bag
Available in blue or red



NITRO SL ROLLATOR (*155) £186.00

Ultra Lightweight Aluminium 5.5kg
Adjustable height handles 32.7-36.6"
Adjustable height backrest
Arthritic friendly, looped cable brakes
Built in stick holder
Removable storage bag
Clip holds frame closed when folded
Max user weight 136kg/21.4st

SATURN ROLLATOR (*199.99) £239.99

Lightweight & modern 7.3kg
8" front soft swivel wheels
Adjustable height handles 34"-41"
Arthritic friendly, looped cable brakes
Built in stick & drink holder
Removable storage bag
Available in black or champagne
Max user weight 22st/136kg



4 WHEEL ROLLATORS (£15*) £18 PER WEEK TO HIRE

TAiMA M-GT (*245.83) £294.99

Ultra lightweight 5.5kg
 As light as a 3-wheeled walker
 Adjustable height handles
 Soft ergonomic handles
 Walking stick holder included
 Large storage bag
 Arthritic friendly, looped cable brake
 Front & side reflectors
 Small & Medium models available
 Max user weight 23kg/150kg



TAiMA XC (*354.16) £424.99
 Pneumatic tyres & XL wheels
 Brilliant for multi surface outdoor use
 Adjustable handle height
 Seat height 61cm
 Total weight 9.14kg
 Max user weight 23st/150kg

EZ Fold-N-Go Walker (*111.99) £134.39

Super lightweight 3.6kg
 2x wheels & 2x easy-glide feet
 No seat & no brakes
 Single handed folding
 Compact & portable
 Storage pouch
 Adjustable hand height 81-97cm
 Available in Blue or Walnut
 Max user weight 28st / 180kg



HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (*20.00) £24.00
HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (*30.00) £36.00

***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**

Try before you buy:
If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.

4 WHEEL ROLLATORS (£15*) £18 PER WEEK TO HIRE



EZ FOLD-N-GO ROLLATOR (*237.99) £285.59

Folds up to 3x smaller than average
Portable & lightweight 6.4kg
Adjustable height handles
Padded handles & locking brakes
2 colours available
Underseat storage bag
Organiser for on the go storage
Max user weight 17.5kg / 113kg

“LETS DREAM” INDOOR ROLLATOR (*319.99) £383.99

Lightweight 6.5kg
Adjustable height single handle
Very narrow: 57cm wide
Full length soft handlebar
Full length one-handed brake
Removable tray
Undertray storage bag
Foldable



HANDY TROLLEY (*85.00) £102.00

Height adjustable handles
Locking brakes on handles
2x Removable trays
Non-folding
Indoor use only
Also available without brakes.

HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (*20.00) £24.00
HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (*30.00) £36.00

***This item is eligible for *VAT relief* if you suffer a chronic medical condition.**

Try before you buy:
If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.

PREMIUM 4 WHEEL ROLLATORS - Not for hire, but can be trialed in the showroom.



TOPRO TROJA 5G (*356.00) £427.20

Back rest included
Adjustable height handles
Available in 4 modern colours/designs
Easy folding, lightweight 7.8kg
Soft grip / ergo grip handles
Huge range of accessories
7 year warranty
3 types of wheels available
2 frame sizes available
Max user weight 23.5st / 150kg.

ROLLZ FLEX 2 (*499.00) £598.80

Stylish, designer rollator, 7.7kg
Large, water resistant, shopping bag
Adjustable angle & height handlebar
Less vibrations from foam-filled tyres
Small turning circle
Leather look seat
Available in 4 colours.
Range of accessories available
Max user weight: 19.5st / 125kg



SALJOL CARBON ROLLATOR (*589.00) £706.80

Ultra Lightweight at 6kg
Integrated drag brake – stops roll-away
All terrain, large, shock absorbing wheels
Cork, natural feel handles
Visible in the dark – over 30 reflectors
Cane holder & removable bag
Available in 3 colours
Large range of accessories available
Seat Heights: 54cm or 62cm
Max user weight: 23.5st / 150kg



SALJOL ALLROUND ROLLATOR (*479.00) £574.80

Very lightweight 6.6kg
Extra large all terrain wheels 25cm
Converts to single handed brake
Upgraded back rest included
Adjustable height natural cork handles
Cane holder comes as standard
Colours: Pearl Grey or Magic Black
Seat Heights: 54cm or 62cm
Max user weight 23.5st / 150kg

