

## WALKING/ZIMMER FRAMES (7.50\*) £9 PER WEEK TO HIRE

Non wheeled from (36.00) £43.20

Wheeled from (40.80) £48.96

Large support base for balance and stability.

Push button height adjustment

PVC hand grips

Easy to wipe clean.

Available with or without 2x front wheels

Folding versions available – from (39.50)



## 4 WHEEL INDOOR ROLLATORS (£15\*) £18 PER WEEK TO HIRE

“LETS GO INDOOR” ROLLATOR (\*264.99) £317.99

Lightweight 6.5kg

Adjustable height single handle

Very narrow: 54.5cm wide

Full length soft handlebar

Full length one-handed brake

Removable tray

Tray colour may vary

Undertray storage bag

Foldable



HANDY INDOOR TROLLEY (\*89.00) £106.80

Height adjustable handles

Locking brakes on handles

2x Removable trays

Non-folding

Indoor use only

Also available without brakes

Safe working load 15kg

Max user weight 20st / 127kg

**HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (\*20.00) £24.00**

**HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (\*30.00) £36.00**

**\*This item is eligible for \*VAT relief\* if you suffer a chronic medical condition.**

**Try before you buy:**

**If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.**

**4 WHEEL ROLLATORS (£15\*) £18 PER WEEK TO HIRE**



**EZ FOLD-N-GO ROLLATOR (\*249.99) £299.99**

- Folds up to 3x smaller than average
- Portable & lightweight 6.4kg
- Adjustable height handles
- Padded handles & locking brakes
- 2 colours available
- Underseat storage bag
- Organiser for on the go storage
- Max user weight 17.5st / 113kg

**EZ Fold-N-Go Walker (\*114.99) £137.99**

- Super lightweight 3.6kg
- 2x wheels & 2x easy-glide feet
- No seat & no brakes
- Single handed folding
- Compact & portable
- Storage pouch
- Adjustable hand height 81-97cm
- Available in blue or walnut
- Max user weight 28st / 180kg



**INDOOR RUTLAND TROLLEY (\*72.50) £87.00**

- Height adjustable 84-99cm
- Two trays
- Non-folding
- No brakes
- Perfect for moving items around the house
- Max user weight 28st / 180kg

**HURRYCANE 2 COMFORT EDITION (\*26.99) £32.39**

- Wide range of walking sticks in stock
- NOT available for hire, all are purchase only



**HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (\*20.00) £24.00**  
**HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (\*30.00) £36.00**

**\*This item is eligible for \*VAT relief\* if you suffer a chronic medical condition.**

**Try before you buy:**

**If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.**