

## **4 WHEEL ROLLATORS (£15\*) £18 PER WEEK TO HIRE**

### **STANDARD ROLLATOR (\*120) £144.00**

Lightweight Aluminium, less than 7kg  
Adjustable height handles 32"-37"  
Adjustable height seat 19"-23"  
Arthritic friendly, looped cable brakes  
Anatomical handgrips  
Padded flip-up seat  
Underseat storage bag  
Available in blue or red  
Max user weight 21<sup>st</sup> / 135kg



### **NITRO SL ROLLATOR (\*155) £186.00**

Ultra Lightweight Aluminium 5.5kg  
Adjustable height handles 32.7-36.6"  
Adjustable height backrest  
Arthritic friendly, looped cable brakes  
Built in stick holder  
Removable storage bag  
Clip holds frame closed when folded  
Available in blue or grey  
Max user weight 136kg/21.4st

### **SATURN ROLLATOR (\*199.99) £239.99**

Lightweight & modern 7.3kg  
8" front soft swivel wheels  
Adjustable height handles 34"-41"  
Arthritic friendly, looped cable brakes  
Built in stick & drink holder  
Removable storage bag  
Available in black or champagne  
Max user weight 22st/136kg



**4 WHEEL ROLLATORS (£15\*) £18 PER WEEK TO HIRE**

**TAiMA M-GT (\*261.66) £313.99**

Ultra lightweight 5.5kg  
As light as a 3-wheeled walker  
Adjustable height handles  
Soft ergonomic handles  
Walking stick holder included  
Large storage bag  
Arthritic friendly, looped cable brakes  
Front & side reflectors  
Small & Medium models available  
Max user weight 23st/150kg



**TAiMA XC (\*429.99) £515.99**

Pneumatic tyres & XL wheels  
Brilliant for multi surface outdoor use  
Adjustable handle height  
Seat height 61cm  
Total weight 9.14kg  
Max user weight 23st/150kg

**TOPRO TROJA 5G (\*356.00) £427.20**

Back rest included  
Adjustable height handles  
Available in 8 modern colours  
Easy folding, lightweight 7.8kg  
Soft grip / ergo grip handles  
Huge range of accessories  
7 year warranty  
2 types of wheels available  
2 frame sizes available  
Max user weight 23.5st / 150kg



**HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (\*20.00) £24.00**

**HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (\*30.00) £36.00**

**\*This item is eligible for \*VAT relief\* if you suffer a chronic medical condition.**

**Try before you buy:**

**If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.**

**4 WHEEL ROLLATORS (£15\*) £18 PER WEEK TO HIRE**

**EZ FOLD-N-GO ROLLATOR (\*249.99) £299.99**

Folds up to 3x smaller than average  
 Portable & lightweight 6.4kg  
 Adjustable height handles  
 Padded handles & locking brakes  
 2 colours available  
 Underseat storage bag  
 Organiser for on the go storage  
 Max user weight 17.5st / 113kg



**EZ Fold-N-Go Walker (\*114.99) £137.99**

Super lightweight 3.6kg  
 2x wheels & 2x easy-glide feet  
 No seat & no brakes  
 Single handed folding  
 Compact & portable  
 Storage pouch  
 Adjustable hand height 81-97cm  
 Available in blue or walnut  
 Max user weight 28st / 180kg

**"LETS GO INDOOR" ROLLATOR (\*264.99) £317.99**

Lightweight 6.5kg  
 Adjustable height single handle  
 Very narrow: 54.5cm wide  
 Full length soft handlebar  
 Full length one-handed brake  
 Removable tray  
 Tray colour may vary  
 Undertray storage bag  
 Foldable



**HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (\*20.00) £24.00**  
**HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (\*30.00) £36.00**

**\*This item is eligible for \*VAT relief\* if you suffer a chronic medical condition.**

**Try before you buy:**

**If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.**

**4 WHEEL ROLLATORS (£15\*) £18 PER WEEK TO HIRE**



**HANDY TROLLEY (\*89.00) £106.80**

Height adjustable handles  
Locking brakes on handles  
2x Removable trays  
Non-folding  
Indoor use only

**RUTLAND TROLLEY (\*72.50) £87.00**

Adjustable height 84-99cm  
Two trays  
Non-folding  
No brakes  
Perfect for moving items around the house  
Max user weight 28st /180kg



**AirGo Fusion 2-in-1 (\*205) £246.00**

Rollator and transport chair  
Easily converts between 2 uses  
Swing-around leg & back rests  
Adjustable height 78.5-92.5cm  
Padded seat & backrest  
Removable storage bag  
67cm width  
Weight 9.5kg  
Max user weight 21st/ 136kg

**HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (\*20.00) £24.00**

**HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (\*30.00) £36.00**

**\*This item is eligible for \*VAT relief\* if you suffer a chronic medical condition.**

**Try before you buy:**

**If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.**

**PREMIUM 4 WHEEL ROLLATORS - Not for hire, but can be trialed in the showroom.**

**ROLLZ FLEX 2 (\*499.00) £598.80**

Stylish, designer rollator, 7.7kg  
Large, water resistant, shopping bag  
Adjustable angle & height handlebar  
Less vibrations from foam-filled tyres  
Small turning circle  
Leather look seat  
Available in 4 colours.  
Range of accessories available  
Max user weight: 19.5st / 125kg



**SALJOL CARBON ROLLATOR (\*589.00) £706.80**

Ultra Lightweight at 6kg  
Integrated drag brake – stops roll-away  
All terrain, large, shock absorbing wheels  
Cork, natural feel handles  
Visible in the dark – over 30 reflectors  
Cane holder & removable bag  
Available in 3 colours  
Large range of accessories available  
Seat Heights: 54cm or 62cm  
Max user weight: 23.5st / 150kg



**SALJOL ALLROUND ROLLATOR (\*479.00) £574.80**

Very lightweight 6.6kg  
Extra large all terrain wheels 25cm  
Converts to single handed brake  
Upgraded back rest included  
Adjustable height natural cork handles  
Cane holder comes as standard  
Colours: Pearl Grey or Magic Black  
Seat Heights: 54cm or 62cm  
Max user weight 23.5st / 150kg



**HIRE ONLY: Zone 1 (up to 20 miles): Delivery & collection charge = (\*20.00) £24.00**  
**HIRE ONLY: Zone 2 (20-40 miles): Delivery & collection charge = (\*30.00) £36.00**

**\*This item is eligible for \*VAT relief\* if you suffer a chronic medical condition.**

**Try before you buy:**  
**If, after the first 2 weeks' hire, you wish to purchase the item, we will deduct the 2 weeks hire & collection charge from the purchase price.**